

NPN uptake assay

Materials

N-Phenyl-1-naphthylamine (NPN), magnesium acetate tetrahydrate were obtained from MACKLIN® (Shanghai, China).

EDTA Na₂ was obtained from Solarbio® (Beijing, China).

Phosphate buffered saline (PBS) was obtained from TransGen Biotech (Beijing, China).

1.5M Tris-HCl was obtained from Beyotime® (Shanghai, China).

Anhydrous Ethanol was obtained from Sinopharm® (Beijing, China).

Hemin was obtained from MedChemExpress® (Shanghai, China).

E. coli ATCC25922 was obtained from hopebio® (Qingdao, China).

Materials needed for DNA origami (see preparation)

Procedures

1. Loading helper strands

- Prepare two microcentrifuge tubes on ice, labeling them as group A and group B. For group A, add M13, S-G4-x, S-x, and S-PAM-cap-x in the indicated amounts. For group B, add M13, S-G4-x, S-x, S-PAM-cap-x, and Apt-cap-x in the indicated amounts. The molar ratio of M13 (5nM) and helper strands (50 nM) is 1:10.

Group A: DNA origami G4 (DO G) (with S-PAM-cap-x)

Reagents	Volume
TAE/Mg ²⁺ /K ⁺ buffer	47μl
S-G4-x	0.25ul × 136 = 34μl
S-x	0.25ul × 50 = 12.5μl
S-PAM-cap-x	0.25ul × 6 = 1.5μl
M13	5μl

S-G4-15	S-G4-66	S-G4-134	S-G4-182	S-133
S-G4-16	S-G4-67	S-G4-135	S-G4-183	S-156
S-G4-17	S-G4-68	S-G4-136	S-G4-184	S-157
S-G4-20	S-G4-69	S-G4-137	S-G4-185	S-180
S-G4-21	S-G4-70	S-G4-138	S-G4-186	S-181
S-G4-22	S-G4-71	S-G4-139	S-G4-187	S-194
S-G4-26	S-G4-74	S-G4-140	S-G4-188	S-198
S-G4-27	S-G4-75	S-G4-141	S-G4-189	S-199
S-G4-28	S-G4-76	S-G4-142	S-G4-190	S-203
S-G4-29	S-G4-77	S-G4-143	S-G4-191	S-204
S-G4-30	S-G4-78	S-G4-146	S-G4-195	S-24

S-G4-31	S-G4-79	S-G4-147	S-G4-196	S-25
S-G4-32	S-G4-80	S-G4-148	S-G4-197	S-48
S-G4-33	S-G4-81	S-G4-149	S-G4-200	S-49
S-G4-34	S-G4-82	S-G4-150	S-G4-201	S-72
S-G4-35	S-G4-83	S-G4-151	S-G4-202	S-73
S-G4-38	S-G4-86	S-G4-152	S-13	S-96
S-G4-39	S-G4-87	S-G4-153	S-14	S-97
S-G4-40	S-G4-88	S-G4-154	S-18	S-120
S-G4-41	S-G4-93	S-G4-155	S-19	S-121
S-G4-42	S-G4-94	S-G4-158	S-23	S-144
S-G4-43	S-G4-95	S-G4-159	S-36	S-145
S-G4-44	S-G4-98	S-G4-160	S-37	S-168
S-G4-45	S-G4-99	S-G4-161	S-60	S-169
S-G4-46	S-G4-100	S-G4-162	S-61	S-192
S-G4-47	S-G4-105	S-G4-163	S-84	S-193
S-G4-50	S-G4-106	S-G4-164	S-85	S-PAM-Cap102
S-G4-51	S-G4-107	S-G4-165	S-89	S-PAM-Cap103
S-G4-52	S-G4-110	S-G4-166	S-90	S-PAM-Cap104
S-G4-53	S-G4-111	S-G4-167	S-91	S-PAM-Cap113
S-G4-54	S-G4-112	S-G4-170	S-92	S-PAM-Cap114
S-G4-55	S-G4-117	S-G4-171	S-101	S-PAM-Cap115
S-G4-56	S-G4-118	S-G4-172	S-108	
S-G4-57	S-G4-119	S-G4-173	S-109	
S-G4-58	S-G4-122	S-G4-174	S-116	
S-G4-59	S-G4-123	S-G4-175	S-125	
S-G4-62	S-G4-124	S-G4-176	S-126	
S-G4-63	S-G4-129	S-G4-177	S-127	
S-G4-64	S-G4-130	S-G4-178	S-128	
S-G4-65	S-G4-131	S-G4-179	S-132	

Group B: DNA origami G4 (DO G) (with S-PAM-cap-x and Apt-cap-x)

Reagents	Volume
TAE/Mg ²⁺ /K ⁺ buffer	44μl
S-G4-x	0.25ul × 136 = 34μl
S-x	0.25ul × 50 = 12.5μl
Apt-cap-x	0.25ul × 12 = 3μl
S-PAM-cap-x	0.25ul × 6 = 1.5μl
M13	5μl

S-G4-15	S-G4-67	S-G4-136	S-G4-185	S-181
S-G4-16	S-G4-68	S-G4-137	S-G4-186	S-194
S-G4-17	S-G4-69	S-G4-138	S-G4-187	S-198
S-G4-20	S-G4-70	S-G4-139	S-G4-188	S-199

S-G4-21	S-G4-71	S-G4-140	S-G4-189	S-203
S-G4-22	S-G4-74	S-G4-141	S-G4-190	S-204
S-G4-26	S-G4-75	S-G4-142	S-G4-191	S-24
S-G4-27	S-G4-76	S-G4-143	S-G4-195	S-25
S-G4-28	S-G4-77	S-G4-146	S-G4-196	S-48
S-G4-29	S-G4-78	S-G4-147	S-G4-197	S-49
S-G4-30	S-G4-79	S-G4-148	S-G4-200	S-72
S-G4-31	S-G4-80	S-G4-149	S-G4-201	S-73
S-G4-32	S-G4-81	S-G4-150	S-G4-202	S-96
S-G4-33	S-G4-82	S-G4-151	S-13	S-97
S-G4-34	S-G4-83	S-G4-152	S-14	S-120
S-G4-35	S-G4-86	S-G4-153	S-18	S-121
S-G4-38	S-G4-87	S-G4-154	S-19	S-144
S-G4-39	S-G4-88	S-G4-155	S-23	S-145
S-G4-40	S-G4-93	S-G4-158	S-36	S-168
S-G4-41	S-G4-94	S-G4-159	S-37	S-169
S-G4-42	S-G4-95	S-G4-160	S-60	S-192
S-G4-43	S-G4-98	S-G4-161	S-61	S-193
S-G4-44	S-G4-99	S-G4-162	S-84	S-PAM-Cap102
S-G4-45	S-G4-100	S-G4-163	S-85	S-PAM-Cap103
S-G4-46	S-G4-105	S-G4-164	S-89	S-PAM-Cap104
S-G4-47	S-G4-106	S-G4-165	S-90	S-PAM-Cap113
S-G4-50	S-G4-107	S-G4-166	S-91	S-PAM-Cap114
S-G4-51	S-G4-110	S-G4-167	S-92	S-PAM-Cap115
S-G4-52	S-G4-111	S-G4-170	S-101	Apt-Cap1
S-G4-53	S-G4-112	S-G4-171	S-108	Apt-Cap2
S-G4-54	S-G4-117	S-G4-172	S-109	Apt-Cap6
S-G4-55	S-G4-118	S-G4-173	S-116	Apt-Cap7
S-G4-56	S-G4-119	S-G4-174	S-125	Apt-Cap11
S-G4-57	S-G4-122	S-G4-175	S-126	Apt-Cap12
S-G4-58	S-G4-123	S-G4-176	S-127	Apt-Cap205
S-G4-59	S-G4-124	S-G4-177	S-128	Apt-Cap206
S-G4-62	S-G4-129	S-G4-178	S-132	Apt-Cap210
S-G4-63	S-G4-130	S-G4-179	S-133	Apt-Cap211
S-G4-64	S-G4-131	S-G4-182	S-156	Apt-Cap215
S-G4-65	S-G4-134	S-G4-183	S-157	Apt-Cap216
S-G4-66	S-G4-135	S-G4-184	S-180	

- b. Mix the above reagents and briefly centrifuge them. The mix is first heated at 95 °C for 10 minutes. Subsequently, annealing is performed by slowly cooling the mixture from 95 °C to 20 °C at a rate of 1 °C/min.
- c. The Amicon® Ultra-0.5 Centrifugal Filter Devices (100kDa) are used to remove redundant helper strands three times. The purified DO G

solution is obtained.

2. Loading functional strands

- a. Add 1.5 μL of the PAM-rich stock solution to group A. Add 1.5 μL of the PAM-rich and aptamer stock solution to group B. Then, restore the volume to 100 μL with TAE/ Mg^{2+} / K^+ buffer.
- b. Mix the above reagents and briefly centrifuge them. The mix is annealed from 45 $^{\circ}\text{C}$ to 25 $^{\circ}\text{C}$ at a rate of 5 min/ $^{\circ}\text{C}$ for six cycles.
- c. The Amicon[®] Ultra-0.5 Centrifugal Filter Devices (100kDa) are used to remove redundant PAM-rich and aptamer three times. The purified $\text{DO}_{\text{PAM}}\text{G}$ and $\text{DO}^{\text{A}}_{\text{PAM}}\text{G}$ solution are obtained.

3. Loading of hemin

- a. Add 1 μL of hemin stock solution to the purified $\text{DO}_{\text{PAM}}\text{G}$ and $\text{DO}^{\text{A}}_{\text{PAM}}\text{G}$ solution obtained in the last step. Then, restore the volume to 100 μL with TAE/ Mg^{2+} / K^+ buffer.
- b. Gently resuspend the solution by pipetting three times using a wide-bore pipette tip.
- c. Briefly centrifuge at $500 \times g$ for 3 seconds to sediment droplets.
- d. Incubate at room temperature for 1 hour.
- e. The Amicon[®] Ultra-0.5 Centrifugal Filter Devices (100 kDa) are used to remove redundant hemin. The purified $\text{DO}_{\text{PAM}}\text{GH}$ and $\text{DO}^{\text{A}}_{\text{PAM}}\text{GH}$ solution are obtained. Then, restore the volume to 100 μL with TAE/ Mg^{2+} / K^+ buffer.

4. Repeat the previous steps for two times.

5. Preparation of NPN stock solution:

- a. Weigh 0.0219g NPN powder using an analytic balanced.
- b. Transfer the weighed powder into 10mL anhydrous ethanol and mix the solution thoroughly. Label the centrifuge tube as NPN-1.
- c. Pipette 1mL NPN-1 into 9mL anhydrous ethanol and mix the solution thoroughly. Label the centrifuge tube as NPN-2.
- d. Pipette 1mL NPN-2 into 19mL PBS and mix the solution thoroughly. Label the centrifuge tube as NPN-3.
- e. Pipette 1mL NPN-2 into 19mL 0.067M Tris-HCl and mix the solution thoroughly. Label the centrifuge tube as NPN-4.
- f. Pipette 1mL NPN-2 into 1mL PBS and mix the solution thoroughly. Label the centrifuge tube as NPN-5.

6. Preparation of EDTA Na_2 solution:

- a. Weigh 0.186g EDTA Na_2 .
- b. Add 100mL ddH₂O to solubilize the weighed EDTA Na_2 .

7. Dilution of bacterial suspension:
 - a. Take 5mL of overnight cultured *E. coli* ATCC25922 suspension and add it into five 1.5mL centrifuge tubes. Centrifuge the five 1.5mL centrifuge tubes at 4000 rpm for 3 minutes.
 - b. Discard the supernatant and resuspended the bacterial sediment in each 1.5mL centrifuge tube with 1mL PBS.
 - c. Pipette 4mL PBS into a 15mL centrifuge tube.
 - d. Add aliquots of bacterial suspension into the 15mL centrifuge tube until the OD600 value approaches 1.5.

8. Preparation of the components of each group:
 - a. Prepare ten 1.5mL centrifuge tubes and label them from 1 to 10.
 - b. For tubes 7-10, add 150uL of bacterial suspension. Centrifuge the four tubes at 12,000 rpm for 3 minutes. Discard the supernatant in each tube. Resuspend the bacterial sediments in tubes 7 and 9 with 28.8μL PBS and resuspend the bacterial sediments in tubes 8 and 10 with 22.8μL PBS.
 - c. Add other components of each group to the corresponding centrifuge tubes according to the table below.

Buffer	1	300μL PBS
	2	240μL PBS 60μL 50μM NPN-3
Buffer + Cells	3	150μL PBS 150μL PBS-diluted bacterial suspension
	4	90μL PBS 60μL 50μM NPN-3 150μL PBS-diluted bacterial suspension
EDTA	5	90μL 0.067M Tris-HCl 60μL 5mM EDTA Na ₂ 150μL PBS-diluted bacterial suspension
	6	30μL 0.2M Tris-HCl 60μL 5mM EDTA Na ₂ 60μL 50μM NPN-4 150μL PBS-diluted bacterial suspension
DO _{PAM} GH	7	267μL 5nM DNA origami 3μL 255μM Hemin 1.2μL 12.5mM H ₂ O ₂
	8	267μL 5nM DNA origami 3μL 255μM Hemin 1.2μL 12.5mM H ₂ O ₂ 6μL 500μM NPN-5
DO ^A _{PAM} GH	9	267μL 5nM aptamer-loaded DNA origami

		3μL 255μM Hemin
		1.2μL 12.5mM H ₂ O ₂
	10	267μL 5nM aptamer-loaded DNA origami
		3μL 255μM Hemin
		1.2μL 12.5mM H ₂ O ₂
		6μL 500μM NPN-5

9. Measurements of fluorescent intensity:

- Transfer the solution in the centrifuge tubes into the black-bottomed 96-well plate. Each group corresponds to three wells on the plate. Incubate the plate at room temperature for 3 minutes.
- Place the black-bottomed 96-well plate into SpectraMax® iD5 Microplate Reader and select “Endpoint” mode under the fluorescence module.
- Set the excitation length to 350nm and the emission length to 420nm.
- Set PMT gain to 500.
- Record the relative fluorescent unit of each well.